


QUICK USER GUIDE

In order to support you when you return home, we have set up remote monitoring of your cardio-respiratory parameters.

As part of this remote monitoring, you have been provided with a **BORA Band BB-100 connected bracelet**, as well as a **BORA Box to run the BORA Connect for Home mobile application**. This document explains how these devices work as well as advice for how to get the best out of them.



Advice for wearing the bracelet

The BORA Band BB-100 should **fit snugly on your wrist** to ensure measurement accuracy. However, you should not feel any discomfort. It must be worn at all times, but should be removed when showering/washing. During this time, we advise you to put it on charge.


Starting up

Press the button for 1 second. The indicator light ● is **static for 5 seconds** and tells you that your BORA Band BB100 is working.

Making sure your BORA Band is working correctly

Press the button for 1 second. The indicator light ● is **static for 5 seconds** and tells you that your BORA Band BB100 is on.





Turning off

To turn off, press the button until the indicator light flashes  for 5 seconds. This means that your BORA Band BB100 is stopping.

Charging

Use the **charger provided by Biosency** to charge your BORA Band. We strongly recommend that you **charge it during those times when you are showering/washing**. Don't forget to put it back on once charged.

Other indicator lights

	Flashing orange light: the battery is low on charge; please charge you BORA Band BB100 using the charger provided.		Blue: The BORA Band is sending or receiving data; please do not turn off during this time.
	Flashing white light: the battery is charging. Once completely charged, this indicator light will switch off.		Fixed red indicator light: your BORA Band BB100 is faulty. Please try switching off and then on again. If there is still a fault, please contact your provider.

BORA Box


BORA Box consists of a smartphone with the “**Bora Connect For Home**” mobile application pre-loaded. This enables it **to transfer data** from your BORA Band to the BORA connect platform.

General usage tips

Leave **the BORA Box connected to its charger** at all times in a room of your house. The telephone must be switched on (the power button can be found on the lower right).

BORA Band communication

The BORA Box automatically connects to your BORA Band when it is switched on and nearby. In the BORA Connect For Home application, you will be able to see the screen opposite, which indicates **the identifier** (or serial number) of your BORA Band, its **battery status**, and **the date of the last data transfer**.



If you do not see this screen, wait for up to 10 minutes. If it does not appear after this time, please go to the “**Transfer problem**” section.

In the event of data transfer problems

If your app does not connect to your BORA Band, please do the following:

1. Make sure your BORA Band is working properly by referring to the “**Checking that the BORA Band is working properly**” section at the beginning of this guide. If this doesn't fix the problem, please continue to the next step.
2. Make sure your BORA Band is close enough to your BORA Box. If this doesn't fix the problem, please contact your service provider.

If the indicator light bothers you during the night

If it is orange (low battery): connect your BORA Band more regularly. For example, we suggest you charge your BORA Band while showering/washing.

If it is blue (data transmission): you can move your BORA Box away from your bedroom to limit data transfers during the night.

For more details, please consult the BORA Band and BORA Connect for Home manuals on the following website:

<https://biosency.com/ifu/>